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STROMA

2.1 Introduction

The Experiment Units STROMA have been used to study the behavior in microgravity of human immune system cells, endothelial cells, mesenchymal stem cells, breast cancer cells and, rodents thyroid cells, muscular cells, osteoclasts and macrophages. The experimental procedures (cell culture, medium changes and cell lysis) on board of the ISS are conducted by the use of specific bioreactors (STROMA) designed by Kayser Italia taking advantage of incubator facilities on-board (KUBIK) or of the ASI BIOKON. The astronaut involvement is limited to the positioning of the bioreactors into the incubator.


2.2 Description of the Experiment Units (EU)

a. Main Feature

The STROMA Experiment Unit is a device capable of performing automatic 2D cell culturing in microgravity. It is equipped with reservoirs for chemicals (culture medium, washing buffer, fixatives, about 1.3 ml of buffers) and a culture chamber allowing growth of adherent cells on a provided support for coverslips. A fluidic path connects reservoirs to the culture chamber allowing fluid injection. The experiment timeline of fully autonomous EU is electronically controlled. At the end of the experiment the STROMA Experiment Unit can be stowed at controlled temperatures, up to -80°C. After stowage and re-entry on Earth, 2D cell cultures can be analysed by microscopy techniques as well as molecular biology-based approaches for genomic, transcriptomic and proteomic studies.

b. Summary of function

Each Experiment Unit is made of a semicrystalline thermoplastic polymer with excellent mechanical and chemical resistance properties, biologically inert. Cross contamination among the chambers is avoided thanks to proper sealing gaskets. The EU itself provides one level of containment (LoC). The experiment is fully autonomous, all the actions are electrically controlled by a pre-defined timeline uploaded on the on board microcontroller. Housekeeping data are recorded during the mission and downloaded at re-entry. On the whole, the actions performed by the fluidic system are achieved by preloaded springs activated by Shape Memory Alloy (SMA) actuators. Such mechanism releases the plungers inward displacing the fluids (Activator or Fixative) contained into the chemicals reservoirs (Activator or Fixative reservoir) toward the Culture Chamber (CC). A manifold channel connects each reservoirs to the CCs so that cells are activated or fixed.

c. Previous mission

The Experiment Units STROMA have been executed on human cells and used in previous space mission, as shown in the following table:

Experiment	Description	Mission /Launch
STROMA	<i>Bone Marrow Stromal Cell Differentiation in Microgravity</i>	2003 with Shuttle STS-107
PITS	<i>Osteoclast cells</i>	2007 with Soyuz-U rocket Foton-M3 capsule
MYO	<i>Muscular cells</i>	2007 with Soyuz-U rocket Foton-M3
SPHINX	<i>Spaceflight of Huvec: an Integrated experiment</i>	2010 with Progress 40P
CYTOSPACE	<i>Microgravity and cells: morphotype and phenotype correlation</i>	2015 with Space X-6
NATO	<i>NANoparticles based countermeasures for Treatment of microgravity induced Osteoporosis</i>	2015 with Space X-6

ENDO	<i>Endothelial cells</i>	2015 Soyuz 44S/Soyuz 42S
SCD	<i>Stem Cell differentiation</i>	2015 with Soyuz 41S
CORM	<i>Retinal Lesion Study</i>	2017 with Space X- 12
NANOROS	<i>Oxidative Stress</i>	2017 with Space X- 12
MYOGRAVITY	<i>Muscle cells</i>	2017 with Space X- 12
NANO	<i>Nano Antioxidants</i>	2019 with Space X-17

ROALD

3.1 Introduction


The Experiment Units ROALD have been used to study the role of apoptosis during human T-lymphocytes activation by mitogens in microgravity and human bone cells.

The experimental procedures (cell culture, culture medium, fixative) board of the ISS are conducted by the use of specific bioreactors (ROALD) designed by Kayser Italia taking advantage of incubator facilities on board (Kubik) or of the ASI BLOKON. The astronaut involvement is limited to the positioning of the bioreactors into the incubator.


3.2 Description of the Experiment Units (EU)

a. Summary of Features

The ROALD Experiment Unit is a device capable of performing automatic cell culture of non-adherent cells (i.e. cells that grow in suspension) in microgravity. It is equipped with reservoirs for chemicals (culture medium, fixatives, about 2.0 ml) and a culture chamber allowing cell growth in suspension. The scientific protocol is led by the KEU-RO electronics following a predefined timeline. At the end of the experiment the ROALD Experiment Unit can be stowed at controlled temperatures, down to -80°C. After stowage and re-entry on Earth, cell cultures can be analyzed with molecular biology-based approaches for genomic, transcriptomic and proteomic studies or cytofluorimetry.

b. Summary of function

Schematic of Insert's Function

Each ROALD Experiment Unit (EU) is made of a semi-crystalline thermoplastic polymer with excellent mechanical and chemical resistance properties, biologically inert. Cross contamination among the fluids chambers are avoided due to proper sealing gaskets. The EU itself provides one *Level of Containment* (LoC) that is increased to two by using KIC-SL containers class. The experiment is fully autonomous; all the actions are electrically controlled by a predefined timeline uploaded into the on-board microcontroller. Housekeeping data are recorded during the mission and downloaded at reentry.

The typical fluidic concept carries out the KEU-RO experimental protocol which relies on three main steps, i.e. cell activation, incubation, fixation. On the whole, the actions performed by the fluidic system are led by preloaded spring actuators activated by the control electronics. Such mechanism pushes the pistons inward displacing the fluids (Activator or Fixative) contained into the chemicals reservoirs (Activator or Fixative reservoir) towards the Culture Chamber (CC). Each CC is provided with a floating piston to allow an expandable volume for fluid injections. Short channels connect independently the reservoirs to the CCs so that cells are activated or fixed (see figures below).

c. Previous mission

The Experiments Units ROALD have been executed and used in previous space mission, as shown in the following table:

<i>Experiment</i>	<i>Description</i>	<i>Mission /Launch</i>
ROALD	<i>ROle of Apoptosis in Lymphocyte Depression</i>	2008 with Soyuz 17S
RESLEM	<i>Role of Endocannabinoid System in human Lymphocytes Exposed to Microgravity</i>	2011 with Soyuz 29S
SERISM	<i>Reprogramming Stem Cells</i>	2017 with Space X-12
BioScience-3	<i>Spaceflight Effects on Vascular Endothelial and Smooth Muscle Cell Processes</i>	2019 with Space X-17